




## Adult breathing systems for transport ventilators




**Critical Care** ▪ Breathing Systems



## Adult breathing systems for transport ventilators

We offer a range of disposable adult transport breathing systems for use with the majority of transport ventilators for inter departmental and inter hospital patient transport or transfer. Set up ready for use, straight from the pack, they use a selection of valves and flow sensors designed and validated for the appropriate ventilators, as shown.

The Oxylog® 3000, 3000 plus and 2000 plus systems also feature an over moulded connector block for the flow sensor lines, ensuring secure connection and correct orientation of the two lines.




System shown is for use with Oxylog® 3000, 3000 plus and 2000 plus.






## HAMILTON-T1




| Code    | Description                                                                           | Box Qty. |
|---------|---------------------------------------------------------------------------------------|----------|
| 2930000 | 30mm Flextube™ transport coaxial breathing system with flow sensor $\geq 1.8\text{m}$ | 12       |
| 2300000 | 22mm Flextube breathing system with flow sensor $\geq 1.6\text{m}$                    | 10       |

## Oxylog® 3000, 3000 plus and 2000 plus



| Code    | Description                                                                                                      | Box Qty. |
|---------|------------------------------------------------------------------------------------------------------------------|----------|
| 2072000 | 22mm Flextube transport breathing system, for use with Oxylog® 3000, 3000 plus and 2000 plus, $\geq 2.0\text{m}$ | 10       |
| 2272000 | 22mm Flextube transport breathing system, for use with Oxylog® 3000, 3000 plus and 2000 plus, $\geq 3.0\text{m}$ | 20       |

## Oxylog® 1000, Pneupac® paraPAC®, ventiPAC™ and transPAC™



| Code    | Description                                                                                                                                                             | Box Qty. |
|---------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| 2070000 | 22mm Flextube basic transport breathing system with non re-breathing valve, for use with Oxylog® 1000, Pneupac® paraPAC®, ventiPAC™ and transPAC™, $\geq 1.2\text{m}$   | 12       |
| 5020002 | 22mm Smoothbore basic transport breathing system with non re-breathing valve, for use with Oxylog® 1000, Pneupac® paraPAC®, ventiPAC™ and transPAC™, $\geq 2.0\text{m}$ | 10       |

illustrated


Oxylog is a registered trademark of Dräger Medical GmbH.

paraPAC and Pneupac are registered trademarks of the Smiths Medical family of companies.


transPAC and ventiPAC are trademarks of the Smiths Medical family of companies.




## Pneupac® paraPAC®, ventiPAC™ and transPAC™



| Code    | Description                                                                                                                                                     | Box Qty. |
|---------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| 2073000 | 22mm Flextube transport breathing system with non re-breathing valve and adjustable PEEP valve, for use with Pneupac® paraPAC®, ventiPAC™ and transPAC™, ≥ 1.2m | 12       |



| Code    | Description                                                                                                                                                          | Box Qty. |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| 2074000 | 22mm Flextube transport breathing system with non re-breathing valve, adjustable PEEP valve and monitoring line with filter, for use with Pneupac® transPAC™, ≥ 1.2m | 15       |



| Code    | Description                                                                                                                                   | Box Qty. |
|---------|-----------------------------------------------------------------------------------------------------------------------------------------------|----------|
| 2075000 | 22mm Flextube transport breathing system with non re-breathing valve and monitoring line with filter, for use with Pneupac® transPAC™, ≥ 1.2m | 15       |

Make an enquiry

Oxylog is a registered trademark of Dräger Medical GmbH.

paraPAC and Pneupac are registered trademarks of the Smiths Medical family of companies.

transPAC and ventiPAC are trademarks of the Smiths Medical family of companies.

IS1.26 • Issue 8 04.22



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016, ISO 14001:2015 and MDSAP

**Please think before you print**  
Save energy and paper.  
If you must print this information sheet please print it double sided.

UK • Ireland • France • Germany • Spain • Portugal • Italy • Benelux • Sweden • Denmark • Lithuania • Russia • Czech Republic  
Turkey • South Africa • China • Japan • Taiwan • Philippines • USA • Canada • Colombia • Australia