


	TITLE: DEVICE MASTER RECORD	
	QUALITY FILE	Ref. no.: QF 04-5
	Revision: 4	Page 1 of 1
	Issue date: 15.07.98	Effective date: 24.10.2024

SUCTION PIECE

REF. NO. **1901 0270 40**



COMPONENT DESCRIPTION

1. SUCTION PIECE

Yankauer Mini, 2 relieving eyes
ID/OD: 2.6/4.0 (CH12), length: 270 mm
Material: Styrolux 684 D
Dye material: PSC, No. 161192, blue

COMPONENT NO.

4.01027040

PACKAGING MATERIAL

Inner bag 40 x 300 mm		2.40300004
Peel pack 81 x 350 mm		5.81350410
Outer bag 260 x 600 mm	100 pcs/bag	2.26060010
Box 386 x 286 x 196 mm	100 pcs/box	6.38628619
Label on the box 85 x 131 mm		7.00085131

Work instruction	QI 03
Packaging instruction	QI 04

Comments: