
	TITLE: DEVICE MASTER RECORD	
	QUALITY FILE	Ref. no.: QF 04-5
	Revision: 2	Page 1 of 1
	Issue date: 24.03.2025	Effective date: 07.11.2025

DuCatheter

REF. NO. **2000 0250 95 US**



COMPONENT DESCRIPTION

COMPONENT NO.

1. SUCTION PIECE

Yankauer "DuCatheter" 4 relieving eyes, 2 x bent
 OD: 9.3 mm, length: 250 mm
 Material: PVC (GC1202D9N2)

0.00025095

PACKAGING MATERIAL

COMPONENT NO.

Pouch 100 x 350 mm
 Outer bag 260 x 600 mm
 Box 318 x 248 x 160 mm
 Label on the box 85 x 131 mm
 Label US on the box 34 x 15 mm

25 pcs/bag
25 pcs/box

5.10035047
 2.26060010
 6.31824816
 7.00085131
 7.00003415

Work instruction
 Packaging instruction

QI 03
 QI 04

Comments: