

DuCatheter by EXTRUDAN®

Oropharyngeal catheter as proposed by Dr. James DuCanto

DuCatheter • 4 RELIEVING EYES • NON STERILE • 2000




Tip Diameter (mm)	(CH)	Length (mm)	Ref-No
6.5 / 9.5	28	250	2000 0250 95 US

DuCatheter • 4 RELIEVING EYES • TUBE PIECE CONNECTOR • NON STERILE • 2004




Tip Diameter (mm)	(CH)	Length (mm)	Ref-No
6.5 / 9.5	28	250	2004 0250 95 US

DuCatheter SUCTION SET • WITH CONNECTING TUBE 90 CM / 180 CM • 4 RELIEVING EYES • NON STERILE • 2000



Tip Diameter (mm)	(CH)	Length (mm)	Tube Diameter (mm)	(CH)	Tube Length (mm)	Ref-No
6.5 / 9.5	28	250	7.0/10.0	30	900	2000 0900 10 ML US
6.5 / 9.5	28	250	7.0/10.0	30	1800	2000 1800 10 ML US

FUNNEL • FUNNEL / FINGERTIP • NON STERILE • 1310 ML FM



Tube Diameter (mm)	(CH)	Length (mm)	Ref-No
7.0/10.0	30	2100	1310 2100 10 ML FM US