

# Blood Collection

## Safety Blood Collection Sets

Safety winged needle with tube and luer adaptor, also available with pre-attached holder. This product provides a closed system for safe venous blood collection.

Once the device is removed, simply slide the safety mechanism over the needle. The click confirms the device is safely locked.

### Features & Benefits

- **Closed system for venous blood collection reduces risk of blood exposure**
- **Tri-cut bevel minimises insertion pain, decreasing patient discomfort**
- **Thin walled for optimal flow rate preserving blood cell integrity**
- **Sharp stainless steel needle for smooth atraumatic insertion**
- **Safety needle reduces the risk of needlestick injury, it easily slides over needle, securely locking it in place**
- **Anti-slip wings for secure handling with gloves**
- **Extension line to provide greater flexibility and comfort for clinician and patient**
- **Luer slip adaptor**
- **Available with pre-attached holder for ease of use**





21G × 3/4" (19mm)

25 pcs / box  
Product Code:

With 20cm Tube, luer adaptor and holder	SBCWN-H-2119-20
With 30cm Tube, luer adaptor and holder	SBCWN-H-2119-30



21G × 3/4" (19mm)

50 pcs / box  
Product Code:

With 20cm Tube and luer adaptor	SBCWN-2119-20
With 30cm Tube and luer adaptor	SBCWN-2119-30



23G × 3/4" (19mm)

25 pcs / box  
Product Code:

With 20cm Tube, luer adaptor and holder	SBCWN-H-2319-20
With 30cm Tube, luer adaptor and holder	SBCWN-H-2319-30



23G × 3/4" (19mm)

50 pcs / box  
Product Code:


With 20cm Tube and luer adaptor	SBCWN-2319-20
With 30cm Tube and luer adaptor	SBCWN-2319-30



25G × 3/4" (19mm)

25 pcs / box  
Product Code:

With 20cm Tube, luer adaptor and holder	SBCWN-H-2519-20
With 30cm Tube, luer adaptor and holder	SBCWN-H-2519-30



25G × 3/4" (19mm)

50 pcs / box  
Product Code:

With 20cm Tube and luer adaptor	SBCWN-2519-20
With 30cm Tube and luer adaptor	SBCWN-2519-30