



Skalpell Engångs Retracting Guard steril

Skalpell Engångs Retracting Guard Steril levereras i förpackningar om 10 st/fp

Sortiment

Skalpell Engångs Steril; 6, 9, 10, 10A, 11, 11P, E/11, Sabre E/11, 12, 12D, 14, 15, 15A, 15C, 15T, Sabre D/15, 16, 18, 19, 20, 21, 22, 22A, 23, Sabre B/23, 24, SGD (Skin Graft Double-edge blade), Stitch Cutters.
Skalpell Engångs Steril med skydd; 10, 10A, 10R, 11, 12, 15, 20, 21, 22, 23 och 24.
Skalpell Engångs Steril Retractable; 10, 11P, 15, 20, 21 samt Stitch Cutter.
Skalpell Engångs Steril Retracting Guard; 10, 11, 11P, 15, 20, 21, 22, 23, 24 samt Stitch Cutter.

Analys knivblad:

ROSTFRITT STÅL

KOL	0.60 - 0.70%
KISEL	0.50% max
MANGAN	1.00% max
KROM	12.00 – 13.50%
SVAVEL	0.025% max
FOSFOR	0.03% max
NICKEL	0.50% max

Material handtag Skalpell Engångs Retracting Guard

Polystyren (High Impact Polystyren)

Material eggskydd Skalpell Engångs Retracting Guard

DuraStar™ Polymer MN611 Natural

Specifikation styckförpackning

Papper: Medicinskt kraftpapper nätlackerad
Film: Polyamid / Polyetylen
Styckförpackningarna är tillverkade inom EU

Förpackning

Kartong: Belagd vitfodrad kartong
Lack: Litho Machine Varnish
Förpackningarna är tillverkade i UK

Styckesförpackning och förpackningar är märkta med aktuell batch och utgångsdatum, förpackningarna är i förekommande fall EAN-märkta.

Produkterna är steriliserade med Etylenoxid av:

Andersen Caledonia, Caledonian House, Phoenix Cres, Bellshill ML4 3NJ USA.



Swann-Morton®
SHEFFIELD ENGLAND



Skalpell Engångs Retracting Guard Steril

<u>Typ</u>	<u>Swann-Morton#</u>	<u>TG Instrument#</u>	<u>Typ</u>	<u>Swann-Morton#</u>	<u>TG Instrument#</u>
10	3601	62330-6SR	20	3606	62337-6SR
11	3603	62331-6SR	21	3607	62339-6SR
11P	3691	62351-6SR	22	3608	62338-6SR
15	3605	62335-6SR	23	3610	62343-6SR
Stitch Cutter	3626	62356-6SR	24	3611	62368-6SR

TG-INSTRUMENT AB

Ingenjörsgatan 9, 252 27 Helsingborg

Vxl 042 12 46 00 - Fax 042 - 12 46 23 - sales@tg-instrument.se - www.tg-instrument.se