

FLOWART® NEEDLE-FREE CHECK VALVE PORT

DESCRIPTION


FlowArt® is a single-use, disposable connector intended for use with a vascular access device for administration of fluids during IV therapy.

FEATURES

- ☑ Compatible with TPN solutions, blood and oncology drugs.
- ☑ Allows flow only towards the patient.
- ☑ ISO 80369-7 compatible.
- ☑ MRI and CT suitable.
- ☑ Shelf Life: 5 Years
- ☑ in compliance with MDD 93/42/EEC. EC Certificate No: 1984-MDD-18-524
- ☑ Available in either bulk / OEM and packed / sterile.
- ☑ Latex-free, PVC-free, DEHP-free

Part	Material
Valve Body	Polycarbonate
Valve Membrane	Silicone

DIAGRAM



PRODUCT DATA

Characteristic	Characteristic
Product Code	A1010CV
Product Name	Needle-Free Check-Valve Port
Priming Volume (mL)	0.30
Flow Rate (mL/min @ 1m water head pressure)	220
Cracking Pressure (mm H ₂ O)	<200
Length (cm)	2.5

LABELLING AND PACKAGING

Each unit is labelled on its blister and packaging container.

Each container is labelled with the following traceability information:

- Quantity,
- Product description,
- Product date,
- Lot number

TDS.003 _Rev02_30.01.2026

Asset® Medikal Tasarım San. Tic. A.Ş. is an ISO 13485 registered company

FlowArt is a registered trademark of Asset Medikal.

Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions.

©2023. All rights reserved.

Asset® Medical, Inc.

5965 Village Way E105-300
San Diego, CA 92130, USA

info@assetmedical.com
www.assetmedical.com

Asset® Medikal Tasarım San. Tic. A.Ş.

İkitelli OSB Mahallesi, 17. Cadde, No:17
Basaksehir, İstanbul - Türkiye

Phone : +90 212 494 27 27
Fax : +90 212 494 34 37