

Product Specifications

Brand : Klinion®
Product : Luer Lock & Luer Slip Infusion Syringes

Product description

Klinion® Luer Infusion Syringes are disposable Luer syringes without needle. They are available in the format of Luer Slip and Luer Lock.

Intended use

Klinion® Luer Infusion Syringes are indicated for injection of fluids for medical purposes into patients of all ages, including neonates, by healthcare professionals. They are compatible with infusion pumps such as Medfusion and Alaris Pump but can also be used for manual injection of fluids for medical purposes.

Classification

Class Is & Im under the Medical Devices Regulations (SI 2002 No. 618, as amended) and
Class Is & Im under Regulation (EU) 2017/745 Medical Device Regulation (MDR).

Packaging & Product Reference

Product Reference	Description	Size (*)	Peel Pouch (pcs)	Inner Carton (pcs)	Outer Carton (pcs)
C0880370	KLINION LUER SLIP SYRINGE 1ML	1 ml	1	100	2400 (24 x 100)
C0880371	KLINION LUER SLIP SYRINGE 3ML	3 ml	1	100	2400 (24 x 100)
C0880372	KLINION LUER SLIP SYRINGE 5ML	5 ml	1	100	2000 (20 x 100)
C0880373	KLINION LUER SLIP SYRINGE 10ML	10 ml	1	100	1200 (12 x 100)
C0880374	KLINION LUER SLIP SYRINGE 20ML	20 ml	1	40	320 (8 x 40)
C0880375	KLINION LUER SLIP SYRINGE 60ML	60 ml	1	30	180 (6 x 30)
C0880380	KLINION LUER LOCK SYRINGE 1ML	1 ml	1	100	2400 (24 x 100)
C0880381	KLINION LUER LOCK SYRINGE 3ML	3 ml	1	100	2400 (24 x 100)
C0880382	KLINION LUER LOCK SYRINGE 5ML	5 ml	1	100	2000 (20 x 100)

Product Specifications

C0880383	KLINION LUER LOCK SYRINGE 10ML	10 ml	1	100	1200 (12 x 100)
C0880384	KLINION LUER LOCK SYRINGE 20ML	20 ml	1	40	320 (8 x 40)
C0880385	KLINION LUER LOCK SYRINGE 60ML	60 ml	1	30	180 (6 x 30)

Product Reference	Description	Size (*)	EAN Peel Pouch	EAN Inner Carton	EAN Outer Carton
C0880370	KLINION LUER SLIP SYRINGE 1ML	1 ml	26971564462642	16971564462706	2697156446270
C0880371	KLINION LUER SLIP SYRINGE 3ML	3 ml	26971564462659	16971564462713	26971564462710
C0880372	KLINION LUER SLIP SYRINGE 5ML	5 ml	26971564462666	16971564462720	26971564462727
C0880373	KLINION LUER SLIP SYRINGE 10ML	10 ml	26971564462673	16971564462737	26971564462734
C0880374	KLINION LUER SLIP SYRINGE 20ML	20 ml	26971564462680	16971564462744	26971564462741
C0880375	KLINION LUER SLIP SYRINGE 60ML	60 ml	26971564462697	16971564462751	26971564462758
C0880380	KLINION LUER LOCK SYRINGE 1ML	1 ml	06971564462648	16971564462645	26971564462642
C0880381	KLINION LUER LOCK SYRINGE 3ML	3 ml	06971564462655	16971564462652	26971564462659
C0880382	KLINION LUER LOCK SYRINGE 5ML	5 ml	06971564462662	16971564462669	26971564462666
C0880383	KLINION LUER LOCK SYRINGE 10ML	10 ml	06971564462679	16971564462676	26971564462673
C0880384	KLINION LUER LOCK SYRINGE 20ML	20 ml	06971564462686	16971564462683	26971564462680
C0880385	KLINION LUER LOCK SYRINGE 60ML	60 ml	06971564462693	16971564462690	26971564462697

*Sizes are approximate

Features

- Graduated
- Accurate dosage injection for pump use.
- Align with Medfusion and Alaris Pump.
- Special dimension design complies with syringe pump injection standard.
- Sterile

Contra-Indications

Klinion® Luer Infusion Syringes are not indicated for multiple-use.

Product Specifications

Precautions

- Single use

Presence of specific components

Component	Present
Latex	No
Phthalates	No
Substances of animal origin	No
Substances of human origin	No
Blood or blood derivatives	No
BPA	No
Heavy metals	No
PVC	No

Sterilization

EO sterilized.

Storage conditions

The products should be stored in a cool, well-ventilated and dry place away from direct sunlight.

Shelf life

Klinion® Luer Infusion Syringes have a shelf-life of 5 years from production date, when stored under specified conditions.

Picture



Product. ref.: C0880385

Product Specifications

Approval	Name	Date	Signature
Assistant Product Manager	Suzanne Tukker		
Associate Category Management Director	Susanne van den Berg		
RA Specialist	Alexandra Maria da Silva Duarte		