



Product Specification

Tourniquet, latex-free

Order number	95.1007.011
Product description	Tourniquet, material: PL, latex-free, 1 piece(s)/bag

Material and Color

Product material	polyester (PL)
Colour of product	red

Product characteristics

Application	venous blood collection
Colour of print/label	white

Purity & Certification

Product category	Medical device (MD)
CE certified	CE supplier
Latex-free	yes
Batched	yes

Packaging

Type of smallest subpackaging	bag
Depth of case	330 mm
Depth of box	80 mm
Width of case	330 mm
Width of box	40 mm
Height of box	40 mm

Packaging

Height of case	220 mm
Case volume	0.024 cbm
Weight of case	6.86 kg
Weight of product	0.0306 kg
Piece (s) / outer case	200
Piece(s) / inner box	1
Piece(s) / pallet	12800
Minimum order qty.	1
EAN case	4260150884806
EAN inner box	4260150884707

This is the current specification for this product. SARSTEDT reserves the right to make changes, in whole or in part, at any time without prior notice.