

Product Specifications

Brand : Klinion®
Product : Fluid Collectors - bags

Product description

The Klinion® Fluid Collectors – bags are a range of plastic bags which include an absorbent pad inside, to collect, transport and dispose of body fluids. The absorbent pad is also available separately. The Klinion® Fluid Collectors are disposable, single-use products.

Intended use

The Klinion® Bedpan Bag, Urinal Bag and Vomit Bag are intended to safely collect, transport, and dispose of body fluids.

The Klinion® Absorbing Pad is intended to quickly absorb liquids and turn them into gels. For use with fluid collectors such as bedpans, urinals, bowls etc.

Classification

Class I under the Medical Devices Regulations (SI 2002 No. 618, as amended) and
Class I under Regulation (EU) 2017/745 Medical Device Regulation (MDR).

Packaging & Product Reference

Product Reference	Description	Size (*)	Inner Carton (pcs)	Outer Carton (pcs)
C150101	Bedpan Bag	62,2 x 45,7 ± 0.5	20	400 (20 x 20)
C150111	Urinal Bag	42,0 x 16,5 ± 0.5	20	400 (20 x 20)
C150121	Vomit Bag	25,0 x 25,0 ± 0.5	20	400 (20 x 20)
C150131	Absorbing Pad	24,0 x 13,0 ± 0.5	30	600 (30 x 20)

Product Specifications

Product Reference	Description	EAN Inner Carton	EAN Outer Carton
C150101	Bedpan Bag	8715343744390	8715343744383
C150111	Urinal Bag	8715343744376	8715343744369
C150121	Vomit Bag	8715343744352	8715343744345
C150131	Absorbing Pad	8715343744338	8715343744321

*Sizes are approximate

Features

- Polyethylene bag with non-woven pad
- Hygienic body waste management
- Strong double drawstring for easy opening and secure closure
- Ready to use

Precautions

- Single use

Specifications

Characteristic	Specification
Material of the bags	LDPE (low density polyethylene)
Seam strength	≥15N/5cm
Pad Absorption	700 ml ± 10% in 90 seconds

Presence of specific components

Component	Present
Latex	No
Phthalates	No
Substances of animal origin	No
Substances of human origin	No
Blood or blood derivatives	No

Product Specifications

Sterilization

Non-sterile.

Storage condition

Klinion® Fluid Collectors should be stored in a dry and well-ventilated room.

Shelf life

The Klinion® Fluid Collectors have a shelf life of 3 years under normal storage conditions.

Picture



Product. ref.:
C150121



Product. ref.:
C150101



Product. ref.:
C150111

Approval	Name	Date	Signature
Associate Category & Manager Facility	Global Product Corina Foucault		
RA Specialist	Alexandra Duarte		