


Product specification

Urine container, max. volume: 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent, with paper label



Product description

Order number	75.9922.745
Product description	Urine container, max. volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, with paper label, label/print: green, 250 piece(s)/bag

Material and Colors

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Size

Length of product	55 mm
Height of product	55 mm
Diameter	44 mm
Length excluding cap	55 mm
Max. volume	70 ml

Product characteristics

Volume of work	50 ml
Base shape	flat base
Label/ Print	with paper label
Colour of print/label	green
Graduation	yes
Closure type	screw cap
Application	collection and storage of urine

Purity & Certification

Product category	In-vitro diagnostic (IVD)
Batched	yes

Packaging

Weight of case	6.6 kg
Piece(s) / pallet	8000
Piece(s) / inner box	250

This is the current specification for this product.
Sarstedt reserves the right to make changes, in full or
in part, at any time without prior notification.

Piece (s) / outer case	500
Width of case	370 mm
Height of case	400 mm
Depth of case	548 mm
Case volume	0.081 cbm
Type of smallest subpackaging	bag
Minimum order qty.	500
EAN inner box	4038917152761
EAN case	4038917027809

This is the current specification for this product.
Sarstedt reserves the right to make changes, in full or
in part, at any time without prior notification.