



Curved Micro Suction

Reference: S446 & S447

Anatomically designed to access harder to reach areas. Perfect for procedures performed through the nose, the smooth rounded tip allows for gentle, atraumatic dissection.

Technical Specification

- 9fg 15cm
- 60° curved tip
- 30° bend
- Single use

Features & Benefits

- Left and right options suitable for both sides of the patients anatomy
- 60° angled tip facilitates use in hard to reach areas
- Available in 'left' (S446) and 'right' (S447) bend options
- Atraumatic rounded tip
- Targeted suction through smaller working spaces
- Safe and robust alternative to difficult to clean re-usable suction
- Lightweight, ergonomic handle with control hole

Order Details

- 25 Pieces per box
- 6 Boxes per packing case



Risk of cross contamination if device re-sterilised.

*Images not to

161215 v4