


X-tendable Tubing Breathing Systems



KEY FEATURES

100% NON PVC

CUSTOM DESIGNS AVAILABLE

Systems for:

ICU, ANAESTHESIA, ENTONOX SYSTEMS

Non-Sterile Clinically Clean Packed

REF CODE	DESCRIPTION	Pack Count Per Carton
22XP1622FF	Basic 1.6M Extended Length	20 units
22XP2022FF	Basic 2M Extended Length	20 units
22DXPG3022FF	2 Limb Breathing System 3M Detachable Y-piece	20 units
22XPG3022FF	2 Limb Breathing System 3M Glued Y-piece	20 units
22DXPG6022FF	2 Limb Breathing System 6M Detachable Y-piece	6 units

Sales Unit 6, 10 & 20 = Transport Cartons

Product Technical Summary

(100% PVC FREE)

Regulatory

[CONTACT US](#)

The product complies with legal requirements for medical devices ISO 11135



REACH regulation:

The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-5, EN ISO 10993-7 and BS EN ISO 10993-10

Compliance with legal requirements: Medical Devices Directive 93/42/EEC, Class IIa

Compliance with product standards: EN ISO 11135, EN ISO 11737-1, BS EN ISO 14971

Quality standards: ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN 1041 and EN ISO 11607-1 , BS EN 15986

Materials: PP/PE Tubing, Connector PP, PE, Styrene

Warnings and Recommendations

Warnings:	The product is single use. If reused performance of product may deteriorate, cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	4 Years
Country of origin:	UK
Product Weight:	300g

Packaging Information

Individual Pack	Material of Pack	Size (L x W x H) mm	Shipping Box	Size (L x W x H) mm
Single Pack Poly Bag	PE	500 x 500mm	Cardboard	L x W x H cm

Barcodes / GTINs

REF	Single Pack	Transport Carton
-----	-------------	------------------