




71.5230 - Pressure monitoring line, blue

Product description

- 1. Pressure monitoring line with a standard bore and blue colour coding to identify venous applications. The tubing has two male luer connectors and is made of PVC with the plasticizer TOTM. Dimensions: 1200 mm, 1.5/2.98 mm (ID/OD), 1.59 ml (tube-volume).



Historical :	
BASIS-UDI-DI :	4046765EXTpvcCLASSICE3



Please note that this document is for information only and neither a specification nor an accompanying document which is relevant to safety. Safety instructions and warnings as well as all other information, which must be observed for the intended use of this product, can be found in the enclosed instructions for use. Subject to alterations without prior notice.

REF	Product Summary	Specific Product Information		UDI
71.5230	Pressure monitoring line, blue		1 100 -	4046765091255 4046765007164 -