

SafeDon™ X-tra

Hygiene System

270 Nitrile

Latex Free • Powder Free • Non Sterile

PRODUCT INFORMATION

Type:	Latex-free, powder-free and non-sterile examination glove
Primary Material:	Nitrile (Acrylonitrile-butadiene)
Design and Feature:	Ambidextrous, straight fingers, textured fingertips, beaded cuff
Colour:	Blue
Powder:	≤ 2 mg/glove, no powder lubricant added
Packing:	200 gloves by weight per box (XL: 180 gloves), 8 boxes per shipping carton
Shelf Life	5 years upon date of manufacture

Dimensions:

Gloves Size	Reorder #	Palm Width (mm)	Min. Length (mm)
Small	SLN270-01	80 ± 10	275
Medium	SLN270-02	95 ± 10	275
Large	SLN270-03	110 ± 10	275
Extra-large	SLN270-04	≥ 110	275

Thickness:

Location of Thickness Measurement	Single Wall (mm)	Typical Value
Finger (at 13 ± 3 mm from the fingertip)	Min. 0.05	0.12
Palm (at center of palm)	Min. 0.05	0.08
Cuff (at 25 ± 5 mm from the cuff edge)	Min. 0.05	0.07

Physical Properties:

Parameters	Before Aging	After Aging
Elongation (%)	Min. 500	Min. 400
Force At Break (N) - specification	Min. 6.0	Min. 6.0
Typical Force At Break (N)	7.0	8.0

Pre-Shipment Inspection:

Inspection Criteria	Inspection Level	AQL
Dimension	S-2	4.0
Physical Properties	S-2	4.0
1000ml Water Leak Test	G-1	1.5

Product Conformance:

Certified as Personal Protective Equipment of Complex Design Category III as per PPE Regulation (EU) 2016/425, type tested in accordance with EN ISO 374-1:2016, EN ISO 374-5:2016 and EN ISO 21420:2020. Medical Device of Class I as per Medical Device Regulation (EU) 2017/745, tested in accordance with EN 455 parts 1-4.

