

# KLINIPLAST fix

## Description

Kliniplast fix range Discontinued, please check Kliniplast fix stretch.





Nonwoven fixing plaster on a roll. Fixation of dressings, catheters and drains. Elastic, soft and comfortable. Easy to cut to size. Nonwoven backing, Polyacrylate adhesive layer.

## Product information

In compliance with: Regulation (EU) 2017/745 Medical Device Regulation (MDR) Class I



## Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3017126	Klinion Kliniplast Fixation Plaster Nonwoven On Roll 10 M X 5 Cm Ref 294340	8 pcs.	192 pcs.			294340
	3017127	Klinion Kliniplast Fixation Plaster Nonwoven On Roll 10 M X 10 Cm Ref 294342	6 pcs.	96 pcs.			294342
	3017128	Klinion Kliniplast Fixation Plaster Nonwoven On Roll 10 M X 15 Cm Ref 294344	4 pcs.	60 pcs.			294344
	3017129	Klinion Kliniplast Fixation Plaster Nonwoven On Roll 10 M X 20 Cm Ref 294346	2 pcs.	52 pcs.			294346