

KLINIFIX elastic

Description






Elastic gauze bandage. White. 58% polyester, 42% viscose. For securing wound dressings, cannulae and padding material. 4m when stretched. Do not use directly on open wounds or damaged skin. Medical Device.

Product information

In compliance with: Regulation (EU) 2017/745 Medical Device Regulation (MDR) Class I



Article overview

| | Article number | Description | Inner carton | Outer carton | EAN (IC) | EAN (OC) | Reference |
|---|----------------|--|--------------|--------------|----------|----------|-----------|
|  | 3052465 | Klinifix Elastic Bandage 4mx4cm 40132227 | 20 pcs. | 960 pcs. | | | 40132227 |
|  | 3052459 | Klinifix Elastic Bandage 4mx6cm 40132228 | 20 pcs. | 840 pcs. | | | 40132228 |
|  | 3052458 | Klinifix Elastic Bandage 4mx8cm 40132229 | 20 pcs. | 600 pcs. | | | 40132229 |
|  | 3052457 | Klinifix Elastic Bandage 4mx10cm 40132230 | 20 pcs. | 480 pcs. | | | 40132230 |
|  | 3052461 | Klinifix Elastic Bandage 4mx12cm 40132231 | 20 pcs. | 360 pcs. | | | 40132231 |