
	Document: PS-032 Product Data Sheet Compression Bag		
	Version: 1.0	Document No: PS-032	Page: 1 of 3
Prepared by: Ketil Thorsen	Approved by: Ketil Thorsen		Date of approval: 22.12.2021




PRODUCT DATA SHEET

Thermosaver Compression bag

Compression bag in water-repellent polyester made for Thermosaver Hypothermic Rescue Bag. The material is the same as the outer layer of the rescue bag itself and it comes in both highviz yellow and grey. The bag has three tightening straps to compress the bag to the maximum. The bag has handles at both ends to get the rescue bag out and in use as quickly as possible. The compression bag also has a shoulder strap that makes it possible to have the rescue bag on your back. The bag can be used several times and can be recycled.

	Document: PS-032 Product Data Sheet Compression Bag		
	Version: 1.0	Document No: PS-032	Page: 2 of 3
Prepared by: Ketil Thorsen	Approved by: Ketil Thorsen		Date of approval: 22.12.2021


	Compression bag
In Use L/W/H cm	70/49/1
Weight g	300
Colour	High Wiz Yellow / Grey
Fabric	100% polyester
Outer skin fabric watercolumn	8000 mm
Shrinkage wash 60deg C	<3%
Biocompatibility lining	Oeko Tex 100
Visibility	ISO 20471

	Document: PS-032 Product Data Sheet Compression Bag		
	Version: 1.0	Document No: PS-032	Page: 3 of 3
Prepared by: Ketil Thorsen	Approved by: Ketil Thorsen		Date of approval: 22.12.2021

Revision History

Revision number	Date	Reason for change	Description of change
1	See published date above	n/a	Initial release

Approval

Approval	Date	Signature
Name: Ketil Thorsen CEO Name: Hiwa Hosseinzadeh OM	22-DEC-2021	KT:  HH: 